


View Brochure :- ECG Machine 303D



Main Features

- ★ Three Channel ECG machine with interpretation
- ★ Compact, portable & light weight
- ★ Simultaneously acquisition of 12 lead ECG Data
- ★ Manual, Auto and Arrhythmia analysis mode
- ★ Rhythm lead for observing abnormal ECG waveform
- ★ Heart rate, automatic measurement and interpretation
- ★ Built in Thermal printer for 80mm
- ★ Real-time 12-ch ECG waveform display
- ★ 16 ECG storage available
- ★ High frequency response upto 160Hz
- ★ USB, RS232 Interface (ECGNET software optional)
- ★ DC/AC power with built-in rechargeable Li-ion battery
- ★ Power supply:110/220V, 50/60Hz



	EC-303D Three Channel ECG with Display
Lead	Standard 12 leads
Acquisition	Simultaneous 12 leads
Resolution	12bit/1000Hz
Operation mode	Manual/Auto/Analysis Modes
Filter	AC,EMG Filter
Drift filter	Anti-Drift System
CMRR	≥ 100dB,with AC Filter
Input circuit	Floating,Protection against defibrillator effect
Input circuit current	≤ 0.1 μ A
Input Impedance	> 50MΩ
Patient current leakage	< 10 μ A
Calibrating voltage	1mV ± 3%
Voltage tolerance	± 500mV
Time constant	> 3.2s
Frequent response	0.05Hz~160Hz
Noise level	≤ 15 μ Vp-p
Threshold	≤ 20 μ V
Paper speed	25,50mm/s(± 3%)
Sensitivity	Auto,2.5,5,10,20,40mm/mV
Recorder	High resolution thermal printer
Safety standard	GB9706.1-2007,IEC Class I,Type CF
Display	320X240
Print paper	80mmX20m
Battery	14.4V(2200mAh)rechargeable Li-ion battery

CONCEPT IMAGING

KD-131, PITAMPURA,
NEAR KOHAT METRO STATION, DELHI-110088
PH: +91-9811239988 / 9312885088

WEBSITE: WWW.CONCEPTIMAGING.IN
EMAIL: SANJEEV_CONCEPT@YAHOO.CO.IN
INFO@CONCEPTIMAGING.IN